

TRANSPORTING THE TRAPEZE

Within a Single Location or Area - The trapeze can be moved by one or two persons, subject to the strength and capability of those person(s). The rear, non-locking wheels are designed to make sliding the trapeze easier. Subject to the area of movement/ distance, the trapeze may be left fully assembled. **ALWAYS** secure the hanging triangle grip before moving the trapeze to avoid injury. To reposition the chain/grip loosen the thumb screw and move the slide ring closer to the mast where the triangle grip may be secured against the mast.

NEVER tip the trapeze on its wheels when moving to avoid possible injury. **ALWAYS** slide the trapeze on the floor to avoid the trapeze falling and causing personal injury, damage to the trapeze or to property.

Disassembling for Storage or to Move to another Location - The trapeze may be fully disassembled for storage or transport. **ALWAYS** pack components with protection materials (e.g., foam pellets) to prevent damage to the trapeze. Follow the assembly instructions in this guide when re-assembling the unit. **ALWAYS** replace any parts that have been damaged during disassembly, transportation or re-assembly.

LIMITED LIFETIME WARRANTY

The steel frame of this ProBasics product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device has been manufactured to exacting standards and carefully inspected prior to shipment. Your ProBasics brand product is warranted to be free of defects in materials and workmanship, has a Limited Lifetime Warranty on the frame and one year warranty from date of purchase on all non-durable parts which are subject to normal wear and need periodic replacement. This warranty provides only for replacement of defective parts and does not cover shipping or labor charges. This warranty is not transferable and only available to end user owners.

This warranty does not cover failure due to owner misuse or negligence, or normal wear and tear. This warranty does not include non-durable components, such as rubber accessories, casters, wheels, chains or grips, which are subject to normal wear and need periodic replacement.

COMPASSHEALTH

Manufactured for:
Compass Health Brands

Distributed by ConvaQuip
Toll-free: 1- 800-637-8436
www.convaquip.com
BariatricTrapeze
Made in China


Free Standing Bariatric Trapeze with Base and Rear Wheels

Model: PBBT

Max User Weight Capacity: 1,000 lb.

*The free-standing trapeze is designed to provide support to the user while changing positions in bed and stability when transferring out of bed. This is an assistive tool only and is **NOT** intended to support the full body weight of the user.*





Product Dimensions

Size	73" h X 43" w X 48.5" d
Leg Width	37.5" Inner/ 43" Outer
Boom Length	36"
Boom Swivel Range	0° - 180°
Product Weight	100 lbs.
Shipping Weight	110 lbs.
Max User Weight Capacity	1,000 lbs.

Map Term	Components	Warranty Definition
Section 1	Left Leg	Frame; wheels are non-durable
Section 2	Right Leg	Frame; wheels are non-durable
Section 3	A Frame	Frame
Section 4	Boom	Frame
A	A Frame Locking Pin/Screw	Non-Durable
B	Slide Ring for Triangle grip Assembly	Non-Durable
C	Slide Ring Locking Thumbscrew	Non-Durable

Tools Required:

Drill, rubber mallet and adjustable wrench

ASSEMBLY INSTRUCTIONS

1. To ensure the legs are positioned correctly, make sure the wheels are located on the inside of the leg bar.
2. Connect sections 1 & 3 by sliding section 3 over 1 and depressing the push pins to secure in place. There should be an audible "click" when the pins are properly engaged. Repeat for section 2.
3. Connect sections 3 & 4 by sliding section 4 down over section 3. Pull the locking pin (A) outward to allow section 4 to lock into position.
4. Slide (B) slide ring over the tube section 4 by depressing push button on section 4 and sliding the into place. Insert the thumb screw (C) through the ring (B) and turn clockwise until tight to secure the ring in place. **NOTE:** thumb screw **MUST** be inserted into the side of B without the receptacle nut to be sure it locks into place correctly.)

IMPORTANT

Prior to patient use, be sure to tighten the trapeze clamps onto the Trapeze base.

OPERATION

Repositioning in Bed or to Assist Getting out of Bed

1. Adjust the height of the triangle so the user may fully and firmly grasp the triangle with both hands without stretching.

OPERATION (CONT'D)

2. Position the triangle grip above the middle of the users' chest by loosening the thumb screw and moving the slide ring. **ALWAYS** tighten the thumb screw before use.
3. Firmly gripping the triangle, the user should pull with both arms evenly, **NEVER** at an angle to help reduce possible injury.

DO NOT use the trapeze to lift or reposition oneself, unless both hands are firmly gripping the triangle and the user is positioned directly below the boom.

WARNING

- **DO NOT** use to support full body weight. This device is intended to assist a user **ONLY** while they are positioned on the bed, **ALWAYS** maintain contact with the bed when using the trapeze. The maximum user weight capacity of this trapeze is 1,000 lbs/ 453 kg. **ALWAYS** observe the weight limits.
- After assembly or transport of the trapeze, **ALWAYS** check that hardware is securely tightened **BEFORE** use.
- **DO NOT** use this device without proper instruction and supervision from a healthcare professional. It is the responsibility of the patient's physician and caregiver to assess the patient and determine which devices and accessories are needed to support and protect the patient, while limiting risk to the patient. Use of fall protection, pressure management, and transfer assistance mechanisms should be considered as appropriate.
- **DO NOT** hang from, stand on, or climb the trapeze stand. **DO NOT** hang objects from the trapeze stand.
- Take care when positioning the trapeze triangle to allow for proper user head clearance.
- **DO NOT** modify or replace parts. If replacement parts are needed, contact an authorized ProBasics dealer.
- **DO NOT** use the trapeze on uneven/rough/sloping ground. Unsuitable surfaces could cause the device to tip.
- **DO NOT** install or adjust the trapeze while the bed is occupied.
- This is not a toy. **DO NOT** let children play on the trapeze.
- **ALWAYS** check tightness of hardware, if the trapeze is disassembled and reassembled for transport. **ALWAYS** use **CAUTION** when moving the trapeze, as moving the unit on uneven surfaces may cause the trapeze to tip.

FAILURE TO USE COMMON SENSE AND HEED THE ABOVE WARNINGS FURTHER INCREASES RISK OF SERIOUS INJURY. USE AT YOUR OWN RISK AND WITH APPROPRIATE AND SERIOUS ATTENTION TO SAFE OPERATION. USE CAUTION.

MAINTENANCE & CARE

EVERY month, inspect all hardware for tightness, fatigue or rust and rubber tips/grips for excessive wear. **ALWAYS** replace the chain, triangle, and/or slide ring/thumb screw if they show signs of excessive wear.

Clean the trapeze triangle at least once a month AND between patients, clean the trapeze stand every 6 months, use mild soap and water, **DO NOT** use abrasive cleaners, allow the trapeze to dry completely before use.