



## BARIATRIC PATIENT SEATING

# BARIATRIC TILT-RECLINE SEATING

The 750-TRC is the choice for tilt and recline positioning of Bariatric patients. Separate manual hand lever controls enable the reclining back and seat tilt to function individually. The seat surface is a metal pan covered with optional two-inch x 28 inch fabric covered high density foam which is removable if needed. Front rigging is Full Width Footboard. Removable Full Length Armrests are standard. Designed for patients weighing up to 750 pounds (340.2 kg).



Back in up position and seat in reverse tilt position.



Reclining back in full upright position and seat in forward tilt position.



Back in full recline position and seat in reverse tilt position.

### Accessories

High Density Foam Seat, Grey Fabric, 2" High Headrest Positioning Wings

### Areas of Use

- Long Term Care
- Respiratory Therapy
- Position Change Therapy
- Patient Room
- Bariatric Care



Push handle and manual control levers. Two levers to control back recline and two levers to control seat tilt.

Back upholstery features adjustable straps to allow upholstery to conform to user

### Specifications:

	Seat Width		Seat Depth		Seat pan to floor		Back Height from seat		Back Angle upright/full incline	forward Seat Tilt	reverse Seat Tilt	Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	degrees	degrees	degrees	lb	Kg	
750-TRC	28	71.1	20	50.8	16.25	41.3	35	88.9	10 - 36	5	22	750	340.2	1/ea

Made in the USA and Distributed by ConvaQuip.

**WARNINGS:** Cancer and Reproductive Harm - [www.P6gWarnings.ca.gov](http://www.P6gWarnings.ca.gov)

Color variations may occur due to monitor or printer settings.

800.637.8426

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.convaquip.com](http://www.convaquip.com)